

baremedical®

Hydrocolloid

Baremedical hydrocolloid dressings aim to provide an optimal wound healing environment with a matrix design that delivers enhanced absorbency, flexibility and overall patient care. Not only do they help prevent residues from leaching into the wound bed, they also keep the contact layer intact, reducing the potential for skin maceration and assists in one-piece removal. The hydrocolloid layer reacts with the cellular material from the wound to form a gel-like covering to help maintain a moist wound environment. The highly absorbent inner layer effectively removes exudate from the wound and facilitates granulation, epithelization and autolysis.

Features & Benefits

- Primary or secondary dressing for wounds with light to moderate exudate
- Low-profile rounded corners and smooth tapered design helps reduce edge roll-up and risk of shear and friction
- Up to 7 day wear time provides undisturbed moist wound environment, minimizes skin maceration and reduces frequency of dressing changes
- Semi-transparent outer layer permits wound site visualization without removing the dressing
- Easy removal system that does not break down or leave residue in the wound



PRODUCT INFORMATION

Sterility	Sterile
Shelf Life	3 years
Latex Status	Latex Free

ORDERING DETAILS

SIZE/STYLE	5 x 5cm	10 x 10cm	15 x 15cm
CH2 CODE	2527914	2527922	2527931
UNIT OF MEASURE	Pack of 20	Pack of 10	Pack of 10